

Product Datasheet

eGFP Antibody (orb1562812)

Description	eGFP Antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Purified recombinant EYFP.tag full length expressed in E.coli.
Target	eGFP
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1.06 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 27 kDa; Observed MW:
Uniprot ID	C8CHS1
Dilution Range	WB: 1/1000-1/10000 IP: 1/20
Expiration Date	12 months from date of receipt.

