

Product Datasheet

Flag Tag (7E1) Mouse mAb Antibody (orb1562815)

Description	Flag Tag (7E1) Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Synthetic Peptide of Flag-Tag
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1.3 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: 1/1000-1/10000 IF: 1/100-1/500 IP: 1/20
Expiration Date	12 months from date of receipt.

