

Product Datasheet

Flag Tag (4A8) Mouse mAb Antibody [Out of stock] (orb1562818)

Description	Flag Tag (4A8) Mouse mAb Antibody [Out of stock]
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide conjugated to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Refer to figures
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.