

Product Datasheet

GFP Antibody (orb1562823)

Description	GFP Antibody
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein of GFP-Tag
Target	GFP
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 27 kDa; Observed MW:
Uniprot ID	Q96319
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.

