

Product Datasheet

Trx Tag (3C1) Mouse mAb Antibody (orb1562844)

Description	Trx Tag (3C1) Mouse mAb Antibody
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant Trx-tag fusion protein expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	2.2 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 12 kDa; Observed MW:
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.

