

Product Datasheet

MAPK14 Antibody (orb1562977)

Description	MAPK14 Antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human P38 MAPK Pprotein fragments expressed in E.coli.
Target	MAPK14
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	2.8 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Uniprot ID	Q16539
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.

