

Product Datasheet

KDM1A Antibody (orb1562979)

Description	KDM1A Antibody
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Purified recombinant human KDM1/LSD1 protein fragments expressed in E.coli.
Target	KDM1A
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	2.5 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 93 kDa; Observed MW: 110 kDa
Uniprot ID	O60341
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Expiration Date	12 months from date of receipt.

