

Product Datasheet

FOXP3 Antibody (orb1563013)

Description	FOXP3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FOXP3.
Target	FOXP3
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Uniprot ID	Q9BZS1
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000
Expiration Date	12 months from date of receipt.

