

Product Datasheet

SLC31A1 Antibody (orb1563021)

Description	SLC31A1 Antibody
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human SLC31A1/CTR1
Target	SLC31A1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Concentration	0.3 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 21 kDa; Observed MW: 26-34 kDa
Uniprot ID	O15431
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.

