

Product Datasheet

KRT18 Antibody (orb1563193)

Description	KRT18 Antibody
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Purified recombinant human Keratin 18 protein fragments expressed in E.coli.
Target	KRT18
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.25 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Uniprot ID	P05783
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200
Expiration Date	12 months from date of receipt.

