

Product Datasheet

GABARAP Antibody (orb1563204)

Description	GABARAP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide.
Target	GABARAP
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.954 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 14 kDa; Observed MW: 14,16 kDa
Uniprot ID	O95166
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.

