

Product Datasheet

HSPA8 Antibody (orb1563223)

Description	HSPA8 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein of Hsc 70.
Target	HSPA8
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Uniprot ID	P11142
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.

