



Product Datasheet

H3C1 Antibody (orb1563253)


Description	H3C1 Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Synthetic Peptide of Histone H3 (Di Methyl Lys27)
Target	H3C1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Uniprot ID	P68431
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of DiMethyl...



Western blot analysis of DiMethyl-Histon...



Immunocytochemistry analysis of DiMethyl-...