

Product Datasheet

BID Antibody (orb1563269)

Description	BID Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	Purified recombinant fragment of human BID expressed in E. Coli.
Target	BID
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Uniprot ID	P55957
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 FC: 1/50-1/100
Expiration Date	12 months from date of receipt.

