

Product Datasheet

S100a10 Antibody (orb1563287)

Description	S100a10 Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Purified recombinant mouse S100A10 protein fragments expressed in E.coli.
Target	S100a10
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	3.4 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Uniprot ID	P08207
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Expiration Date	12 months from date of receipt.

