

Product Datasheet

EIF4EBP2 Antibody (orb1563292)

Description	EIF4EBP2 Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human 4E-BP2.
Target	EIF4EBP2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	5.66 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 13 kDa; Observed MW: 15-20 kDa
Uniprot ID	Q13542
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.

