

Product Datasheet

CCND2 Antibody (orb1563388)

Description	CCND2 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human Cyclin D2 protein fragments expressed in E.coli.
Target	CCND2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1.9 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	Calculated MW: 33 kDa; Observed MW: 38 kDa
Uniprot ID	P30279
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.

