

Product Datasheet

Acetyl-Lysine Rabbit pAb Antibody (orb1563427)

Description	Acetyl-Lysine Rabbit pAb Antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Acetylated BSA
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.

