

Product Datasheet

AAPK1/AAPK2 Antibody (orb1563453)

Description	AAPK1/AAPK2 Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human AMPK alpha around the phosphorylation site of Thr172.
Target	AAPK1/AAPK2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Uniprot ID	Q13131 , P54646
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 ELISA: 1/10000
Expiration Date	12 months from date of receipt.

