

Product Datasheet

A1BG Antibody (orb1563489)

Description	A1BG Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from A1BG.
Target	A1BG
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 54 kDa; Observed MW: 28 kDa
Uniprot ID	P04217
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000
Expiration Date	12 months from date of receipt.

