

Product Datasheet

MAP2 Antibody (orb1563493)

Description	MAP2 Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	Synthetic Peptide of MAP2
Target	MAP2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P11137
Dilution Range	IHC: 1/50-1/100 IF: 1/50-1/200
Expiration Date	12 months from date of receipt.

