

Product Datasheet

CHEK1 (Phospho-Ser296) antibody (orb156350)

Description	CHEK1 (Phospho-Ser296) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CHEK1 isoform b around the phosphorylation site of Ser296 (IQ(p-S)NL)
Target	Phospho-CHEK1 (Ser296)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	O14757
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.