

Product Datasheet

FTL Antibody (orb1563531)

Description	FTL Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant full length of human ferritin light chain protein expressed in E.coli.
Target	FTL
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	4.7 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P02792
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.

