

Product Datasheet

c-Kit(Phospho-Tyr823) antibody (orb156379)

Description	c-Kit(Phospho-Tyr823) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human c-Kit around the phosphorylation site of Tyr823 (SN(p-Y)VV)
Target	phospho-c-Kit (Tyr823)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	145 kDa
Uniprot ID	P10721
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.