


Product Datasheet

CLSTN1 antibody (orb156407)

Description	CLSTN1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Calsyntenin 1 (501-600/981 aa)
Target	CLSTN1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89/107 kDa
Uniprot ID	O94985
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of (mo)nephrocy...