



Product Datasheet

CNG channel beta 1 antibody (orb156412)

Description	CNG channel beta 1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human truncated GARP (851-950/1251 aa)
Target	CNGB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140 kDa
Uniprot ID	Q14028
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse brain tissue



Immunohistochemical staining of Mouse brain tissue