


Product Datasheet

Complement C4 antibody (orb156427)

Description	Complement C4 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human Complement C4 C4a anaphylatoxin or complement C4 alpha chain (701-800/1744 aa)
Target	Complement C4 / C4a
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8/84/190 kDa
Uniprot ID	P0C0L4
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...