

Product Datasheet

Coproporphyrinogen III Oxidase antibody (orb156438)

Description	Coproporphyrinogen III Oxidase antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CPOX/Coproporphyrinogen Oxidase (361-454/454 aa)
Target	Coproporphyrinogen III Oxidase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	P36551
Expiration Date	12 months from date of receipt.