



Product Datasheet

COQ10B antibody (orb156439)

Description	COQ10B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human COQ10B (38-110/238 aa)
Target	COQ10B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	Q9H8M1
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:HepG2(Human) Cell;HL60(Huma...



Western blot:Lung (Mouse);Primary:Anti-C...