

Product Datasheet

C-REL (Phospho-Ser503) antibody (orb156451)

Description	C-REL (Phospho-Ser503) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human C-REL around the phosphorylation site of Ser503 (TS(p-S)DS)
Target	phospho-C-REL (Ser503)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69 kDa
Uniprot ID	Q04864
Dilution Range	WB=1:500-2000, ICC=1:25, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.