


Product Datasheet

CTLA4 antibody (orb156486)

Description	CTLA4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CTLA-4/CD152 (75-170/223 aa)
Target	CTLA4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	P16410
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Western blot analysis of k562 lysates us...