

Product Datasheet

C-TYQEDRKFYKIRD antibody (orb156493)

Description	C-TYQEDRKFYKIRD antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from Zea mays of viviparous-8 (251-350/713 aa)
Target	viviparous8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78 kDa
Uniprot ID	A0A1D6L709
Expiration Date	12 months from date of receipt.