


Product Datasheet

DAGLA antibody (orb156533)

Description	DAGLA antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human DAGLA
Target	DAGLA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	115 kDa
Uniprot ID	Q9Y4D2
Dilution Range	IHC-P: 1:100-500, IHC-F: 1:100-500, ICC: 1:100-500, IF: 1:100-500, ELISA: 1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...