



Product Datasheet

Tlr2 Antibody (orb1565339)


Description	Cooperates with LY96 to mediate the innate immune response to bacterial lipoproteins and other micro...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	Recombinant protein of mouse TLR2
Target	Tlr2
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	0.3 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Affinity Purified
MW	Calculated MW: 89 kDa
Uniprot ID	Q9QUN7
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Expiration Date	12 months from date of receipt.



Western blot analysis of TLR2 in Raw264....



Immunocytochemistry analysis of TLR2 (gr...



Immunohistochemistry analysis of paraffi...