

Product Datasheet

CTNND1 Antibody (orb1565376)

Description	CTNND1 Antibody
Species/Host	Rabbit
Reactivity	Hamster, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Tyr228 of human delta 1 Catenin
Target	CTNND1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Concentration	0.3 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 108 kDa; Observed MW: 95,100 kDa
Uniprot ID	O60716
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Expiration Date	12 months from date of receipt.

