

Product Datasheet

DAXX (Phospho-Ser668) antibody (orb156541)

Description	DAXX (Phospho-Ser668) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human DAXX around the phosphorylation site of Ser668 (LP(p-S)PP)
Target	phospho-DAXX (Ser668)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81 kDa
Uniprot ID	Q9UER7
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.