

Product Datasheet

SUMO1 Antibody (orb1565559)

Description	SUMO1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, FC, IHC-P, IP, WB
Immunogen	A synthetic peptide of human Sumo 1
Target	SUMO1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 12 kDa; Observed MW: 80,12 kDa
Uniprot ID	P63165
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20 FC: 1/50-1/100 ChIP: 1/20
Expiration Date	12 months from date of receipt.

