

Product Datasheet

**delta 1 Catenin (Phospho-Thr310) antibody
(orb156568)**

Description

delta 1 Catenin (Phospho-Thr310)

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human delta 1 Catenin around the phosphorylation site of Thr31 (TG(p-T)PS)
Target	phospho-delta 1 Catenin (Thr310)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	108 kDa
Uniprot ID	O60716
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.