



Product Datasheet

DNA Polymerase gamma antibody (orb156612)

Description	DNA Polymerase gamma antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human POLG (1101-1239/1239 aa)
Target	POLG
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140 kDa
Uniprot ID	P54098
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human placenta



Immunohistochemical staining of human colon