

Product Datasheet

DNASE1L1 antibody (orb156621)

Description	DNASE1L1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human DNASE1L1 (201-302/302 aa)
Target	DNASE1L1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	P49184
Dilution Range	Flow-Cyt=2ug/Test, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.