

Product Datasheet

LPIN1 Antibody [Out of stock] (orb1566211)

Description	LPIN1 Antibody [Out of stock]
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC-Fr, IHC-P, IP, WB
Immunogen	A synthetic peptide of human Lipin 1
Target	LPIN1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 99 kDa; Observed MW: 130 kDa
Uniprot ID	Q14693
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100
Expiration Date	12 months from date of receipt.