

Product Datasheet

ADAM10 Antibody [Out of stock] (orb1566217)

Description	ADAM10 Antibody [Out of stock]
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human ADAM10
Target	ADAM10
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 84 kDa; Observed MW: 90,68 kDa
Uniprot ID	O14672
Dilution Range	WB: 1/500-1/1000 IP: 1/20
Expiration Date	12 months from date of receipt.