


Product Datasheet

DOCK2 antibody (orb156626)

Description	DOCK2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human DOCK2 (451-550/1830 aa)
Target	DOCK2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	212 kDa
Uniprot ID	Q92608
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...