

Product Datasheet

ACTB Antibody (HRP) (orb1566573)

Description	ACTB Antibody (HRP)
Species/Host	Mouse
Reactivity	Algae, Drosophila, Goat, Hamster, Human, Monkey, Mouse, Plant, Rabbit, Rat, Saccharomyces, Yeast, Zebrafish
Conjugation	HRP
Tested Applications	WB
Immunogen	Purified recombinant human beta-actin protein fragments expressed in E.coli.
Target	ACTB
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	2.4 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Uniprot ID	P60709
Dilution Range	WB: 1/5000-1/10000
Expiration Date	12 months from date of receipt.

