

Product Datasheet

COX4I1 Antibody (orb1566591)

Description	COX4I1 Antibody
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IHC-P, IP, WB
Immunogen	A synthetic peptide corresponding to carboxyl terminal residues of human COX IV
Target	COX4I1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	4.6 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 20 kDa; Observed MW: 17 kDa
Uniprot ID	P13073
Dilution Range	WB: 1/1000 IHC: 1/50-1/200 IF: 1/100-1/200 IP: 1/20 FC: 1/100
Expiration Date	12 months from date of receipt.

