

Product Datasheet

Oxfendazole (orb1566850)

Description	Oxfendazole
Reactivity	Other
Tested Applications	ELISA, WB
Preservatives	0.2M PBS, pH7.2, 10% glycerol
Concentration	1mg/ml
Storage	This product is stable for up to twelve months from date of receipt at 4°C, Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	$C_{15}H_{13}N_3O_3S$
Source	This compound was modified in house and tested by animal assays and ELISA.
CAS Number	53716-50-0
Expiration Date	12 months from date of receipt.