

Product Datasheet

EEF2 (Phospho-Thr56 + Thr58) antibody (orb156701)

Description	EEF2 (Phospho-Thr56 + Thr58) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EEF2 around the phosphorylation site of Thr56+Thr58 (RF(p-T)D(p-T)RK)
Target	Phospho-EEF2 (Thr56 + Thr58)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	95 kDa
Uniprot ID	P13639
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.