

Product Datasheet

EEF2k (Phospho-Ser359) antibody (orb156702)

Description	EEF2k (Phospho-Ser359) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EEF2k around the phosphorylation site of Ser359 (CG(p-S)P)
Target	Phospho-EEF2k (Ser359)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82 kDa
Uniprot ID	O00418
Expiration Date	12 months from date of receipt.