


Product Datasheet

Eph receptor B1 (Phospho-Tyr928) antibody (orb156749)

Description	Eph receptor B1 (Phospho-Tyr928) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EphB1 around the phosphorylation site of Tyr928 (VQ(p-Y)RD)
Target	phospho-Eph receptor B1 (Tyr928)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	106 kDa
Uniprot ID	P54762
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...