

Product Datasheet

EPO Receptor (Phospho-Tyr485) antibody (orb156762)

Description

EPO Receptor (Phospho-Tyr485)

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human EPO Receptor around the phosphorylation site of Tyr485 (GP(p-Y)SN)
Target	phospho-EPO Receptor (Tyr485)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	P19235
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.